



T-01 Forces in a 2-bar Truss

(1 **TS-01** Truss member Set1; 1 **SL-01** & **SL-02** String Lines with end clasps)

Aims: This experiment is designed to allow students to observe measure and compare the static displacements and associated force equilibrium conditions at the common connection node of a 2-bar truss (both members of equal initial length and stiffness) to loading applied at that node (orientation of load to be selected by the instructor from a Table).

Learning Outcomes:

After performing this experiment, students will be able to:

- (i) Use the Spreadsheet supplied to compare the predicted forces in each member of a 2-bar truss induced by the applied load using the method of joints with the experimentally obtained values of these forces.
- (ii) Identify the influence orientation of the applied load relative to the orientation of each member in the two-bar truss would have on the induced internal member forces.
- (iii) Reinforce their understanding of the equilibrium of co-planar forces in the context of the equilibrium of forces at the pin-joint of a 2-bar truss.
- (iv) Appreciate the role played by a member's characteristic stiffness in the stiffness method in the evaluation of internal actions from the individual member elongation/contraction.

Equipment/Resources Required:

- (i) **TM-00** (Pixi with window frame in “landscape” configuration with transparent film & pens; Set of chrome/stainless steel bearing balls (weight forces) and load buckets; Digital Scales for weight force evaluation)
- (ii) A digital camera – Better than 8Mega Pixel preferred (Hi-Res phone cameras are suitable)
- (iii) **xyRectify** photogrammetric software on a suitable Notebook, Laptop or PC – if desired.

